



**Reliable Dental Laboratory Dallas**

# **Cerec Connect Step by Step Instructions**

**“Meeting All Your Dental Lab Needs”**

10610 Metric Dr. Suite #195  
Dallas, TX 75243  
(972) 272.5511 Office  
(214) 503.8686 Fax  
[www.rdlldental.com](http://www.rdlldental.com)

## **GETTING STARTED**

**Verify each of these steps before using CEREC Connect:**

- Check the Internet connection on the Acquisition unit
- Download the software onto the Acquisition unit
- Register with CEREC Connect

## **SENDING CASES WITH CEREC CONNECT**

**Step by Step Instructions:**

1. From your desktop screen click the CEREC Connect icon
2. Click add new patient
3. Click ADD NEW CASE
4. Select tooth # and restoration type
5. Proceed to acquisition phase
6. Take images of the appropriate arch in the proper folder
7. Orient the quadrants so that their angle matches the buccal bite
8. Grab the buccal bite and drag it toward either quadrant, matching up images where they overlap. It will flip to orient itself if necessary.
9. The matched images will show a mottled pattern where they overlap. Next, drag the merged images onto the opposing arch.
10. When the three images are matched, the buccal bite will appear mottled on both the preparation and antagonist and click NEXT.
11. Click NEXT to move to the die trimming step.
12. Position the lines for sectioning the model, then click NEXT.
13. Trace the margin and click NEXT.
14. Enter patient and restoration information and click NEXT STEP.
15. Select **“Reliable Dental Laboratory Dallas”** from the drop-down menu. Select Return Date and Return Time and click NEXT.
16. Enter any additional information you wish the lab to have such as “call doctor to discuss case.” Click ADD TO CART.
17. The case data will begin to upload to CEREC Connect. The time will vary depending on the upload capacity of your Internet connection and the file size.
18. Make sure that the correct case is showing in the drop-down menu, then click NEXT STEP.
19. Verify final case information and click SUBMIT. After clicking SUBMIT you have completed your CEREC Connect case submission.

