



# Reliable Dental Laboratory

**“Meeting All Your Fixed, Removable & Dental Implant Needs”**

## COSMETIC & IMPLANT CHECKLIST

*To help serve you better, we have put together a checklist needed  
to complete a cosmetic or extensive case*

- 1. Model of pre-op
- 2. Model of diagnostic wax-up
- 3. Photos of pre-op teeth with shade guide from different angles
- 4. Photos of preps with stump shade tabs
- 5. Bite registration - Right, Left, Anterior/segments of the provisional left in place to maintain vertical while taking the bites in a tripod fashion
  - a) A bite taken with the upper provisional in place against the lower preparations
  - b) A bite taken with the lower provisional in place against the upper preparations
- 6. Stick bite, facial plane or Kois bow (horizontal plane parallel to interpupillary line)
- 7. Photos of temporaries, full face including eyes
- 8. Photos of temporaries, eye to chin relaxed lip position
- 9. Model of prototype (provisionals) and additional comments
- 10. Indicate the desired length of final restorations
- 11. Detailed lab Rx expressing goals of patient
- 12. *With all Implant cases* indicate type and size of each implant to be restored
- 13. *With Implant cases* indicate your preference to the type of part to be utilized
- 14. For multiple model work or multiple dies you *must* indicate which dies are to be utilized

**Email photos to: [photo@rldental.com](mailto:photo@rldental.com)**

### **Additional photos that are helpful:**

- Full Face
- Profile
- Rest Position
- Relaxed Smile
- E Smile
- Retracted View

**“Working Together To Create Smiles For A Lifetime”**

10610 Metric Dr. #191 • Dallas, TX 75243 • Office 972.272.5511 Or Toll Free 855.808.5511

[www.reliabledentallab.com](http://www.reliabledentallab.com)